



Opto[®] Plasma B Discovery Sample Prep Kit, Mouse

PRODUCT INSERT

Purpose

The Opto[®] Plasma B Discovery Sample Prep Kit, Mouse (750-02050) consists of reagents required to screen mouse plasma B cells with the Opto[®] B Discovery workflow on the Beacon[®] optofluidic systems.

Kit Contents

Table 1 – Box 1 (750-08108) components to be stored at 4°C

Component	Unit Quantity	Volume/ Unit	Uses/Unit	Storage Conditions	Part Number
Basal Media	1	185 mL	1	4°C	750-02050
Media Additive A2	1	220 µL	1		
Wetting Additive ^a	2	1.5 mL	1		

a. This reagent is also provided as part of the Wetting Additive Kit (750-08091).

Table 2 – Box 2 (750-08109) components to be stored at -20°C

Component	Unit Quantity	Volume/ Unit	Uses/Unit	Storage Conditions	Part Number
DNA Clean Up ^a	1	20 µL	1	-20°C	750-02050
Loading Reagent ^b	1	1 mL	1		
Media Additive A1	1	30 mL	1		
Media Additive P1 (M)	1	220 µL	1		
Media Additive P2 (M)	1	220 µL	1		

a. This reagent can be purchased separately (520-70000).

b. This reagent is also provided as part of the Loading Reagent Kit (750-08029).

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

Instructions for Use

Instructions for use and handling can be found in the following documents on the Support Portal:

- *Opto® B Discovery Workflow User Manual* (MAN-08351)
- *Plasma B Cell Isolation for Mouse Tissues Protocol* (MAN-000048)

Safety

Use standard laboratory safety protocols and always wear appropriate personal protective equipment (PPE). Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, email support@brukercellularanalysis.com.

Legal

The Opto Plasma B Discovery Sample Prep Kit, Mouse product insert and its contents are proprietary to Bruker Corporation. This product insert is intended to be used, and should only be used, in connection with the use of the Beacon optofluidic systems. No other use is permitted without the prior written consent by an authorized officer of Bruker. No right, title, or interest in any Bruker intellectual property rights are conveyed or granted, expressly or by implication, by this product insert. Review the documents in [Instructions for Use](#) PRIOR to using the Opto Plasma B Discovery Sample Prep Kit, Mouse.

Bruker reserves the right to change the information in this product insert without notice.

Bruker Cellular Analysis
brukercellularanalysis.com
+1 (510) 858-2855
sales@brukercellularanalysis.com

Technical Support
support.brukercellularanalysis.com
+1 (888) 254-5595 ext. 4
support@brukercellularanalysis.com

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

Bruker Cellular Analysis, Beacon, Opto, and the Bruker logo are trademarks of Bruker Corporation. All other marks are the property of their respective owners.