



SpotLight™ Human Lambda Reagent

PRODUCT INSERT

Purpose

The SpotLight™ Human Lambda Reagent (520-08130) consists of the assay reagent required to run a single 4-chip Opto® CLD workflow with up to three assays per chip (twelve total) on the Beacon® optofluidic system, or two 2-chip Opto CLD workflows with up to three assays per chip on the Beacon Select™ optofluidic system. This reagent binds the Lambda light chain of human antibodies and is labeled with a red fluorescent dye.

Kit Contents

| Component | Unit Quantity | Volume/ Unit | Uses/ Unit | Storage Conditions | Part Number |
|-------------------------------------|---------------|--------------|------------|--------------------|-------------|
| SpotLight Human Lambda Reagent (8X) | 1 | 400 µL | 1-2 | 4°C | 520-08130 |

Instructions for Use

Instructions for use and handling can be found in the following documents on the Support Portal:

- *Beacon® Optofluidic Systems User Manual (MAN-000003)*
- *Opto® CLD 2.0 Workflow User Manual (MAN-08096)*

Safety

Use standard laboratory safety protocols and always wear appropriate personal protective equipment (PPE). Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, email support@brukercellularanalysis.com.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

Legal

The SpotLight Human Lambda Reagent product insert and its contents are proprietary to Bruker Corporation. This product insert is intended to be used, and should only be used, in connection with the use of the Beacon optofluidic systems. No other use is permitted without the prior written consent by an authorized officer of Bruker. No right, title, or interest in any Bruker intellectual property rights are conveyed or granted, expressly or by implication, by this product insert. Review the documents in [Instructions for Use](#) PRIOR to using the SpotLight Human Lambda Reagent.

Bruker reserves the right to change the information in this product insert without notice.

Bruker Cellular Analysis
brukercellularanalysis.com
+1 (510) 858-2855
sales@brukercellularanalysis.com

Technical Support
support.brukercellularanalysis.com
+1 (888) 254-5595 ext. 4
support@brukercellularanalysis.com

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

Bruker Cellular Analysis, Beacon, Beacon Select, Opto, OptoSelect, SpotLight, and the Bruker logo are trademarks of Bruker Corporation. All other marks are the property of their respective owners.